

MAXPAR® DN3 Antibody Labeling Kits

The following series of MAXPAR® DN3 Antibody Labeling Kits is available. Each kit contains polymer, tag isotope, buffers, washes and instructions. The 1X kit is sufficient to conjugate up to 100 µg of customer-supplied antibody; the 10X kit is sufficient to conjugate up to 1000 µg of customer-supplied antibody. Each kit includes the enriched stable isotope tag indicated, and the tag isotopes are interchangeable. These kits are designed for labeling reduced disulfides in the Fc fragment of the antibody, through a bismaleimide linker. Please contact the company for other needs.

Product Description	1X Catalog #	10X Catalog #
MAXPAR® DN3 Antibody Labeling Kit, 139La (*)	201139C	201139D
MAXPAR® DN3 Antibody Labeling Kit, 141Pr	201141C	201141D
MAXPAR® DN3 Antibody Labeling Kit, 142Nd	201142C	201142D
MAXPAR® DN3 Antibody Labeling Kit, 143Nd	201143C	201143D
MAXPAR® DN3 Antibody Labeling Kit, 144Nd	201144C	201144D
MAXPAR® DN3 Antibody Labeling Kit, 145Nd	201145C	201145D
MAXPAR® DN3 Antibody Labeling Kit, 146Nd	201146C	201146D
MAXPAR® DN3 Antibody Labeling Kit, 147Sm	201147C	201147D
MAXPAR® DN3 Antibody Labeling Kit, 148Nd	201148C	201148D
MAXPAR® DN3 Antibody Labeling Kit, 149Sm	201149C	201149D
MAXPAR® DN3 Antibody Labeling Kit, 150Nd	201150C	201150D
MAXPAR® DN3 Antibody Labeling Kit, 151Eu	201151C	201151D
MAXPAR® DN3 Antibody Labeling Kit, 152Sm	201152C	201152D
MAXPAR® DN3 Antibody Labeling Kit, 153Eu	201153C	201153D
MAXPAR® DN3 Antibody Labeling Kit, 154Sm	201154C	201154D
MAXPAR® DN3 Antibody Labeling Kit, 156Gd	201156C	201156D
MAXPAR® DN3 Antibody Labeling Kit, 158Gd	201158C	201158D
MAXPAR® DN3 Antibody Labeling Kit, 159Tb	201159C	201159D
MAXPAR® DN3 Antibody Labeling Kit, 160Gd	201160C	201160D
MAXPAR® DN3 Antibody Labeling Kit, 162Dy	201162C	201162D
MAXPAR® DN3 Antibody Labeling Kit, 164Dy	201164C	201164D
MAXPAR® DN3 Antibody Labeling Kit, 165Ho	201165C	201165D
MAXPAR® DN3 Antibody Labeling Kit, 166Er	201166C	201166D
MAXPAR® DN3 Antibody Labeling Kit, 167Er	201167C	201167D
MAXPAR® DN3 Antibody Labeling Kit, 168Er	201168C	201168D
MAXPAR® DN3 Antibody Labeling Kit, 169Tm	201169C	201169D
MAXPAR® DN3 Antibody Labeling Kit, 170Er	201170C	201170D
MAXPAR® DN3 Antibody Labeling Kit, 171Yb	201171C	201171D
MAXPAR® DN3 Antibody Labeling Kit, 172Yb	201172C	201172D
MAXPAR® DN3 Antibody Labeling Kit, 174Yb	201174C	201174D
MAXPAR® DN3 Antibody Labeling Kit, 175Lu (*)	201175C	201175D
MAXPAR® DN3 Antibody Labeling Kit, 176Yb	201176C	201176D

(*) note: In, La and Lu are natural-abundance labels.

In is 95.7% ¹¹⁵In and 4.3% ¹¹³In.

La is 99.9% ¹³⁹La and 0.09% ¹³⁸La.

Lu is 97.4% ¹⁷⁵Lu and 2.59% ¹⁷⁶Lu.

