

MAXPAR® X8 Antibody Labeling Kits

The following series of MAXPAR® X8 Antibody Labeling Kits is available. Each kit contains polymer, tag isotope, buffers, washes and instructions. The 1X kit is sufficient to conjugate up to 100 µg of customer-supplied antibody; the 10X kit is sufficient to conjugate up to 1000 µg of customer-supplied antibody. Each kit includes the enriched stable isotope tag indicated, and the tag isotopes are interchangeable. These kits are designed for labeling reduced disulfides in the Fc fragment of the antibody, through a bismaleimide linker. Please contact the company for other needs.

Product Description	1X Catalog #	10X Catalog #
MAXPAR® X8 Antibody Labeling Kit, 139La (*)	201139A	201139B
MAXPAR® X8 Antibody Labeling Kit, 141Pr	201141A	201141B
MAXPAR® X8 Antibody Labeling Kit, 142Nd	201142A	201142B
MAXPAR® X8 Antibody Labeling Kit, 143Nd	201143A	201143B
MAXPAR® X8 Antibody Labeling Kit, 144Nd	201144A	201144B
MAXPAR® X8 Antibody Labeling Kit, 145Nd	201145A	201145B
MAXPAR® X8 Antibody Labeling Kit, 146Nd	201146A	201146B
MAXPAR® X8 Antibody Labeling Kit, 147Sm	201147A	201147B
MAXPAR® X8 Antibody Labeling Kit, 148Nd	201148A	201148B
MAXPAR® X8 Antibody Labeling Kit, 149Sm	201149A	201149B
MAXPAR® X8 Antibody Labeling Kit, 150Nd	201150A	201150B
MAXPAR® X8 Antibody Labeling Kit, 151Eu	201151A	201151B
MAXPAR® X8 Antibody Labeling Kit, 152Sm	201152A	201152B
MAXPAR® X8 Antibody Labeling Kit, 153Eu	201153A	201153B
MAXPAR® X8 Antibody Labeling Kit, 154Sm	201154A	201154B
MAXPAR® X8 Antibody Labeling Kit, 156Gd	201156A	201156B
MAXPAR® X8 Antibody Labeling Kit, 158Gd	201158A	201158B
MAXPAR® X8 Antibody Labeling Kit, 159Tb	201159A	201159B
MAXPAR® X8 Antibody Labeling Kit, 160Gd	201160A	201160B
MAXPAR® X8 Antibody Labeling Kit, 162Dy	201162A	201162B
MAXPAR® X8 Antibody Labeling Kit, 164Dy	201164A	201164B
MAXPAR® X8 Antibody Labeling Kit, 165Ho	201165A	201165B
MAXPAR® X8 Antibody Labeling Kit, 166Er	201166A	201166B
MAXPAR® X8 Antibody Labeling Kit, 167Er	201167A	201167B
MAXPAR® X8 Antibody Labeling Kit, 168Er	201168A	201168B
MAXPAR® X8 Antibody Labeling Kit, 169Tm	201169A	201169B
MAXPAR® X8 Antibody Labeling Kit, 170Er	201170A	201170B
MAXPAR® X8 Antibody Labeling Kit, 171Yb	201171A	201171B
MAXPAR® X8 Antibody Labeling Kit, 172Yb	201172A	201172B
MAXPAR® X8 Antibody Labeling Kit, 174Yb	201174A	201174B
MAXPAR® X8 Antibody Labeling Kit, 175Lu (*)	201175A	201175B
MAXPAR® X8 Antibody Labeling Kit, 176Yb	201176A	201176B

(*) note: In, La and Lu are natural-abundance labels.

In is 95.7% ¹¹⁵In and 4.3% ¹¹³In.

La is 99.9% ¹³⁹La and 0.09% ¹³⁸La.

Lu is 97.4% ¹⁷⁵Lu and 2.59% ¹⁷⁶Lu.

